


1. Part Number: 2210A-1000

2. Photo:



3. Product Information:

Description	Cohesive Slot Sampler 1000mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Outer Diameter of Outer Tube:	25mm (plus 10mm approx. for scrapers)
Inner Diameter of Inner Tube:	17mm
Overall Length:	1025mm
Slot Length:	120mm
Slot Width:	15mm
Number of Slots	4
Slot Positions:	



External Polish:	< 0.8 Microns Ra
Internal Polish:	None
Nominal Weight of Sampler:	2.0Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Oct-20