


**1. Part Number:** 2210A-1500

**2. Photo:**



**3. Product Information:**

Description	Cohesive Slot Sampler 1500mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Outer Diameter of Outer Tube:	25mm (plus 10mm approx. for scrapers)
Inner Diameter of Inner Tube:	17mm
Overall Length:	1550mm
Slot Length:	120mm
Slot Width:	15mm
Number of Slots	5
Slot Positions:	



External Polish:	< 0.8 Microns Ra
Internal Polish:	None
Nominal Weight of Sampler:	2.7Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 7-Oct-20