

1. Part Numbers: A264-10

2. Photograph:



3. Product Information:

Description	Stainless Bucket 10L
Nominal Volume:	10L
Brimful Volume:	10.5L
Material of Construction:	304 stainless steel
Method of Construction:	Crevice free interior
Surface Finish:	Brush Finish - better than 1.0 microns Ra
Max. Diameter:	300mm
Bucket Height:	258mm
Graduations:	Yes (every 1L)
Skirt Band:	Yes
Nominal Weight:	1.3 kg
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jun-20