

Product Specification

Specification Number

8126H-1000A

Revision

04

Part No. 8126H-1000A

 Description Wide Mouth HDPE Bottle 1000ml Amber
 Materials of Construction:

 Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)
 Cap PP (conforms to FDA 177.1520 & USP Class VI)

 Manufacturing Environment: Controlled environment within factory
 Nominal Volume: 1000ml
 Brimful Volume: 1075ml
 Nominal Height: 201mm
 Diameter: 92mm
 Internal Neck Diameter: 52mm
 Weight: 102g
 Pack Size: 6 bottles per pack
 Method of Sterilisation: N/A (these bottles are non-sterile)
 BSE/TSE Free: Yes
 Supplied with lid screwed on: Yes

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18