

Product Specification

Specification Number

8126H-1000

Revision

04

Part No. 8126H-1000

Description Wide Mouth HDPE Bottle 1000ml

Materials of Construction:

Bottle	HDPE	(conforms to FDA 177.1520 & USP Class VI)
Cap	PP	(conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment: Controlled environment within factory

Nominal Volume: 1000ml

Brimful Volume: 1075ml

Nominal Height: 201mm

Diameter: 92mm

Internal Neck Diameter: 52mm

Weight: 100g

Pack Size: 6 bottles per pack


Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes

Supplied with lid screwed on: Yes

Specification Approved By:

Signature:



Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18