

Product Specification

Specification Number

8126H-125A

Revision

03

Part No. 8126H-125A

Description Wide Mouth HDPE Bottle 125ml Amber
 Materials of Construction:

Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)
 Cap PP (conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment: Controlled environment within factory

Nominal Volume: 125ml

Brimful Volume: 140ml

Nominal Height: 99mm

Diameter: 51mm

Internal Neck Diameter: 28mm

Weight: 27g

Pack Size: 12 bottles per pack

Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes

Supplied with lid screwed on: Yes

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18