

Product Specification

Specification Number

8126H-250

Revision

04

Part No. 8126H-250

Description Wide Mouth HDPE Bottle 250ml

Materials of Construction: Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)
 Cap PP (conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment: Controlled environment within factory

Nominal Volume: 250ml

Brimful Volume: 290ml

Nominal Height: 131mm

Diameter: 62mm

Internal Neck Diameter: 32mm

Weight: 37g

Pack Size: 12 bottles per pack

Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes

Supplied with lid screwed on: Yes

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18