

Product Specification

Specification Number

8126H-30A

Revision

02

Part No. 8126H-30A

Description Wide Mouth HDPE Bottle 30ml Amber

Materials of Construction:

Bottle	HDPE	(conforms to FDA 177.1520 & USP Class VI)
Cap	PP	(conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment: Controlled environment within factory

Nominal Volume: 30ml

Nominal Height: 63mm

Diameter: 37mm

Internal Neck Diameter: 21mm

Weight: 13g

Pack Size: 12 bottles per pack

Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes

Supplied with lid screwed on: Yes

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18