

Product Specification

Specification Number

8126H-500

Revision

04

Part No. 8126H-500

Description Wide Mouth HDPE Bottle 500ml

Materials of Construction: Bottle HDPE (conforms to FDA 177.1520 & USP Class VI)
Cap PP (conforms to FDA 177.1520 & USP Class VI)



Manufacturing Environment: Controlled environment within factory

Nominal Volume: 500ml

Brimful Volume: 555ml

Nominal Height: 169mm

Diameter: 73mm

Internal Neck Diameter: 43mm

Weight: 68g

Pack Size: 12 bottles per pack

Method of Sterilisation: N/A (these bottles are non-sterile)

BSE/TSE Free: Yes

Supplied with lid screwed on: Yes

Specification Approved By:

Signature:

Name: Jonathan Fenner

Position: Director

Date: 9-Jan-18