


	<h1>Product Specification</h1>	Specification Number	8131A-50
		Revision	01

Part No. 8131A-50  
Description Sample Containers 50ml



Nominal Volume 50ml  
Brimful Volume 70ml  
Colour: Translucent with red lid  
Label Area: 50 x 25mm  
Graduation Marks: None  
Material of Construction: Body Polypropylene (PP)  
Lid High Density Polyethylene (HDPE)  
Overall Height (with lid on): 59mm  
Internal Neck Diameter: 43mm  
Body Height: 55mm  
Nominal Weight: 12g  
Individually Bagged? No  
Pre-sterilised? No  
Number of Containers per Box: 384  
Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 14-Feb-14
----------------------------	--	--