

Part No. 8888P-250S

Description: PharmaBottle® 250ml

Use: Single use, disposable bottle

Material of Construction: Polypropylene (PP) – opaque

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011



Method of Construction: Single piece, injection moulded

Height: 88mm

Lid Diameter: 76mm

Nominal Weight: 57g

Assembly & Packing Environment: Class 100,000 Medical Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of containers per Box: 25

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaBottle are registered trade marks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 04-Dec-12
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	