

	<h1>Product Specification</h1>	Specification Number	A736-1000
		Revision	01

Part No. A736-1000

Description Beaker (316L stainless steel)

Nominal Volume: 1000ml

Material of Construction: 316L stainless steel



Method of Construction: Crevice free, all welds ground and polished


Surface Finish: Better than 0.5 microns Ra

Height: 145mm

Internal Diameter: 100mm

Overall Diameter: 110mm

Nominal Weight: 480g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 25-May-16
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	