


	<b>Product Specification</b>	Specification Number	A736-2000
		Revision	01

Part No. A736-2000  
Description Beaker (316L stainless steel)  
Nominal Volume: 2000ml  
Material of Construction: 316L stainless steel



Method of Construction: Crevice free, all welds ground and polished  
Surface Finish: Better than 0.5 microns Ra  
Height: 155mm  
Internal Diameter: 140mm  
Overall Diameter: 155mm  
Nominal Weight: 780g

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 25-May-16
Page 1 of 1	<b>Sampling Systems Ltd.</b> 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	