


**1. Part Number:** 2220B-1500

**2. Photo:**



### 3. Product Information:

Description	Cohesive Pocket Sampler 1500mm
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Overall Length:	1550mm
Outer Diameter of Outer Tube:	25mm (plus 10mm approx. for scrapers)
Inner Diameter of Inner Tube:	17mm
Pocket Length:	135mm
Pocket Width:	11mm
Pocket Volume:	20ml
Number of Pockets	5
Pocket Positions:	



External Finish:	< 0.8 Microns Ra
Internal Pocket Finish:	Fine machined finish
Nominal Weight of Sampler:	4.7 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Oct-20