

1. Part Number: 8370P-200

2. Photograph:



3. Product Information:

Description:	Liquid Funnel 200mm
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Max. Diameter:	232mm
ID at Top of Funnel:	200mm
Overall Height:	180mm
OD at Bottom of Spout:	19mm
ID at Bottom of Spout:	15mm
OD at Top of Spout:	33mm
Nominal Weight:	92g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 28-Feb-23