

1. Part Number: A105-150

2. Photograph:



3. Product Information:

Description	PTFE Funnel 150mm
Material of Construction:	PTFE
Construction:	Single piece, crevice free
Max. Outer Diameter:	158 mm
Outer Diameter of Spout:	18 mm
Inner Diameter of Spout:	11 mm
Length of Spout:	78 mm
Overall Height:	195 mm
Nominal Weight:	300 g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Nov-19