

## **Specification Sheet**

Specification Number: A105-50
Revision: 02
Date Issued: 8-Nov-19

**1. Part Number:** A105-50

## 2. Photograph:



## 3. Product Information:

Description PTFE Funnel 50mm

Material of Construction: PTFE

Construction: Single piece, crevice free

Max. Outer Diameter: 53 mm

Outer Diameter of Spout: 10 mm

Inner Diameter of Spout: 5 mm

Length of Spout: 50 mm

Overall Height: 82 mm

Nominal Weight: 23 g

Quantity per Pack: 1

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora
Position: Quality Manager

Date: 8-Nov-19

