

1. Part Number: A105-70

2. Photograph:



3. Product Information:

Description	PTFE Funnel 70mm
Material of Construction:	PTFE
Construction:	Single piece, crevice free
Max. Outer Diameter:	80 mm
Outer Diameter of Spout:	12 mm
Inner Diameter of Spout:	7 mm
Length of Spout:	60 mm
Overall Height:	110 mm
Nominal Weight:	52 g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Nov-19