


 <b>SAMPLING SYSTEMS</b>	<h1>Sampler Specification</h1>	Specification Number	8051H-1100
		Revision	02

Part No. 8051H-1100

Description Disposable Jumbo Pipette – 1100mm

Material of Construction:

Tip (handle)	HDPE (virgin, FDA acceptable CFR 177.1520 & EU 10/2011)
Outer Tube	HDPE (virgin, FDA acceptable CFR 177.1520 & EU 10/2011)
Tip (product end)	HDPE (virgin, FDA acceptable CFR 177.1520 & EU 10/2011)



Method of Construction: Ultrasonic welding

Nominal Length: 1100mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: Handle End (Blue) = 1.5mm  
Product End (Natural) = 3mm

Nominal Weight of Sampler: 120g

Assembly & Packing Environment: Class 100,000 Medical Cleanroom


Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 50

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 29-Nov-12
Page 1 of 1	<b>Sampling Systems Ltd.</b> 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	