


	<h1 style="margin: 0;">Sampler Specification</h1>	Specification Number	8051H-1100
		Revision	04

Part No. 8051H-1100

Description Disposable Jumbo Pipette – 1100mm

Material of Construction: High Density Polyethylene (HDPE) - virgin

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004



Method of Construction: Ultrasonic welding

Nominal Length: 1100mm

Diameter of Outer Tube: 21mm

Diameter of hole in tip: Handle End (Blue) = 1.5mm
Product End (Natural) = 3mm

Nominal Weight of Sampler: 120g

Assembly & Packing Environment: Class 100,000 Cleanroom


Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 37

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 27-Jul-18
----------------------------	--	--