


	<h1>Product Specification</h1>	Specification Number	8131A-100
		Revision	01

Part No. 8131A-100
Description Sample Containers 100ml



Nominal Volume 100ml
Brimful Volume 145ml
Colour: Translucent with red lid
Label Area: N/A
Graduation Marks: Every 20ml
Material of Construction: Body Polypropylene (PP)
Lid High Density Polyethylene (HDPE)
Overall Height (with lid on): 76mm
Internal Neck Diameter: 56mm
Body Height: 72mm
Nominal Weight: 14g
Individually Bagged? No
Pre-sterilised? No
Number of Containers per Box: 280
Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 14-Feb-14
----------------------------	--	--