

	<h1>Product Specification</h1>	Specification Number	8888P-500
		Revision	01

Part No. 8888P-500

Description: PharmaBottle® 500ml

Use: Single use, disposable bottle

Material of Construction: Polypropylene (PP) – opaque

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011



Method of Construction: Single piece, injection moulded

Height: 106mm

Lid Diameter: 94mm

Nominal Weight: 83g

Assembly & Packing Environment: Class 100,000 Medical Cleanroom

Individually Bagged? Yes (heat sealed PE bag)


Method of Sterilisation: N/A

Number of containers per Box: 25

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare and PharmaBottle are registered Trademarks belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 19-Dec-12
----------------------------	----------------------------------------------------------------------------------------------------	----------------------------------------------------------------