


	<h1>Product Specification</h1>	Specification Number	A416-5
		Revision	02

Part No. A416-5
Description Stainless Container 5L

Materials of Construction:

- Body: 316 stainless steel
- Lid: 316 stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)




Method of Construction: All welds ground and polished
Crevice free interior

Surface Finish: Better than 0.5 microns Ra

Dimensions
A = 119mm
B = 98mm
C = 316mm
D = 164mm

Weight: 2.7 Kg (body + lid + gasket + clamp)

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 4-Jul-14
----------------------------	----------------------------------------------------------------------------------------------------	---------------------------------------------------------------