

	<h1>Product Specification</h1>	Specification Number	A836-1000
		Revision	02

Part No. A836-1000


Description Stainless Bottle 1000 mL

Nominal Volume: 1000 mL

Brimful Volume 1200 mL

Material of Construction:

Body:	316 stainless steel
Lid:	316 stainless steel
Clamp:	304 stainless steel
Gasket:	Silicone (conforms to FDA CFR 177.2600)



Method of Construction: Crevice free interior, welds ground & polished


Surface Finish: Better than 0.5 microns Ra

Dimensions:

A = 64 mm
B = 48 mm
C = 178 mm
D = 102 mm

Nominal Weight (body, lid, clamp & gasket): 980 g

Note: The neck of the container is fitted with a 2" diameter ferrule which conforms to BS4825-3

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 17-May-16
Page 1 of 1	Sampling Systems Ltd. 4-5 Forge Mills Park, Station Road, Coleshill, Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Fax: +44(0)1675 466994 Email: info@sampling.com Web: www.sampling.com	